




# PORTABLE RAMPS

Rubber Ramps  
Dec-Pac Ramps  
Bifold Ramps  
Telescopic Ramps  
Multifolding Ramps

**Intelligent Mobility Solutions**

# PORTABLE RAMPS

## Telescopic Ramps



The Telescopic Ramp allows wheelchair users to mount stairs easily and safely. It is ideal for wheelchair users requiring the assistance of any attendant. The telescopic ramp extends and shortens in length to suit the required height. The non-slip tread provides greater traction and the underside of ramp and is covered with a non-slip rubber surface material. The ramp also has safety catch locks to secure the ramp stays together when not in use. The nylon carrying strap is to assist transportation.

### Features

- Suitable for wheelchairs and scooters
- Maximum load of 180kg on both ramps
- Lightweight robust construction
- Gritted surface for increased grip
- Extremely easy to use and transport

## Dec-Pac Ramps



Dec-Pac Ramps are available in 4 model, in a wide range of sizes:

**PERSONAL** - Portable, lightweight ramp that allows easier access in and around the home or workplace

**MULTIPURPOSE** - Longer and wider, designed to cope with a wider range of obstacles making it more suitable for commercial and travel situations

**SENIOR** - Flexible range of ramps to cope with obstacles from 22.5-87cm in height

**EDGE BARRIER LIMITER (EBL)** - Added security of edge barriers, particularly suited to providing instant access to buses, trains and buildings where a permanent ramp is impractical.

### Features

- Quickly folded into a single compact unit
- Made from a unique combination of uni-directional glass and carbon fibre materials
- Multi-paneled design held together with a patented hinge system connected by webbing
- Tested to AS/NZ Standards, up to 300kg
- No assembly required

**SIZES AVAILABLE FROM 70 to 350CM!**

# PORTABLE RAMPS

## Rubber Ramps



Rubber Ramps have been developed for ease of access through doorways and entrances into shower recesses. These ramps have superior stability and are slip-resistant, reducing the risk of trips and falls.

### Features

- Suitable for shower recesses
- **Stable** - will not slip or curl reducing the risk of trips and falls
- **Weather resistant** - perfect for outdoor use and impervious to liquids
- Tough and long lasting
- **Lightweight** - easy to handle and install
- Made from 100% recycled rubber

**MADE TO YOUR ORDER!**



## Multifolding Ramps

Multifolding Ramps offer the length required for wheelchairs and scooters to access minivans, full-size vans, steps, and raised landing with ease. Max. user capacity of 272kg.

### Features

- Separates into two pieces for easy carrying
- Folds into Half in width and half in length carries like a suitcase
- Patented centre hinge joint eliminates a pinch point and adds strength
- Durable welded aluminium fabrication
- Anti slip traction surface
- Lightweight easy to handle and set up

## Bifold Ramps

Bifold Ramps are ideal for wheelchair users requiring the assistance of an attendant. They are made from aluminium and max. user capacity is 200kg. It is easy to fold, storage and transport them.



### Features

- Ideal for wheelchairs requiring the assistance of an attendant
- Folds for storage and transportation
- Carry handle to assist transportation
- Aluminium construction



**Active Mobility**

**Active Mobility Systems**

Unit 8, 110-120 Silverwater Rd., SILVERWATER NSW 2128

**P:** 02 9649 2111    **F:** 02 9649 8506

**W:** [www.activemobility.com.au](http://www.activemobility.com.au)    **E:** [sales@activemobility.com.au](mailto:sales@activemobility.com.au)

