



Active Mobility



ALIVIO

Masterpiece

LIFT CHAIR COLLECTION

A COMBINATION OF STYLE AND SUPPORT WITH COMFORT

www.activemobility.com.au
sales@activemobility.com.au
02 9649 2111

PICASSO

Dual Motor Lift Chair



Description and Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy-to-use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower head position



FOOTREST
EXTENSION

Specifications

Overall Depth	1003mm
Overall Weight	813mm
Overall Height	1080mm
Seat Depth	508mm
Seat Width	521mm
Seat Height	495mm
Reclined Length	1830mm
Product Weight	74kg
Max User Weight	158kg

LEONARDO

Dual Motor Lift Chair



Description and Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy-to-use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower head position



LAY-FLAT
RECLINER

Specifications

Overall Depth	1067mm
Overall Weight	902mm
Overall Height	1092mm
Seat Depth	521mm
Seat Width	571mm
Seat Height	483mm
Reclined Length	1830mm
Product Weight	76kg
Max User Weight	158kg

MICHELANGELO

Dual Motor Lift Chair



Description and Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy-to-use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower head position



LUMBAR SUPPORT TO ACCOMMODATE
SITTING POSTURE

Specifications

Overall Depth	1080mm
Overall Weight	890mm
Overall Height	1067mm
Seat Depth	546mm
Seat Width	559mm
Seat Height	470mm
Reclined Length	1830mm
Product Weight	70kg
Max User Weight	158kg

DONATELLO

Dual Motor Lift Chair



Description and Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy-to-use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower head position

Specifications

Overall Depth	1016mm
Overall Weight	890mm
Overall Height	1092mm
Seat Depth	520mm
Seat Width	520mm
Seat Height	508mm
Reclined Length	1830mm
Product Weight	74kg
Max User Weight	158kg



STARNDARD OR
HIGH SUPPORT
BACKREST



LATERAL SUPPORT
BACKREST

DONATELLO PETITE

Dual Motor Lift Chair



Description and Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy-to-use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower head position

Specifications

Overall Depth	978mm
Overall Weight	851mm
Overall Height	1041mm
Seat Depth	457mm
Seat Width	457mm
Seat Height	457mm
Reclined Length	1830mm
Product Weight	70kg
Max User Weight	158kg



STARNDARD
SUPPORT BACKREST



LATERAL SUPPORT
BACKREST

Four independent motors for a wide range of position combinations and support



Power Headrest

The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.



Power Lumbar

Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.



Power Backrest & Footrest

With a the dual lift system, it allows you to individually adjust the backrest and the footrest for the position that best fits the clinical needs of the client.



Easy to use hand control:

- Soft grip buttons
- 6 pre-defined positions
- Choose "full recline" or "full stand" with the touch of a button



The **two-stage lift mechanism** has a more vertical lift with a slight tilt compared to the forward tilt you find on generic chairs.

General Conditions that the Power Headrest and Lumbar adjustment can assist with:

Kyphosis · Lordosis · Arthritis



Customer Feedback:

I came over to your shop yesterday to try out the latest chair to hit Australia and wow what a chair! There is nothing in your shop the comes close and I've tried them all. Having a fused neck, arthritis and just had two lower back surgeries I am needing support while sitting. I could never get a chair to support my lower back and neck until now because of the electric adjustable lumbar support and electric adjustable head rest. Gone are the days of trying to use pillows to get the support which could never be achieved. Even the way the chair lifts you has improved as you don't feel like you are going to land on your face.

Therapist Feedback:

I feel this would have a good application for kyphosed clients and offer adjustability not in the typical lift chairs ranges on offer. The contoured seat base is also a good feature, to assist with lower limb and pelvic positioning, acting as a lateral thigh guide to a degree. Although often client with a higher or more complex seating needs may not be the best candidate for the standard lift chairs ranges, these will be welcome additions to provide low to moderate seating support and correction which I think will be beneficial.



Improving your life everyday™

www.activemobility.com.au
sales@activemobility.com.au
02 9649 2111