

C1000 Transactive bariatric ceiling hoist

Bariatric clients often exceed a hoist's safe working load and this can create an unsafe working environment. The C1000 hoist is the strongest hoist of its class with a SWL of 454Kg, making it possible to lift or transfer patients of all sizes with respect and dignity.

With an option of four styles of carry bars, all featuring the Quick Release System, consideration can now be given to the inevitable health related issues that bariatric patients have and they can be lifted according to their needs. The C1000 comes with a digital display with built-in diagnostics and lift counter for user maintenance and compliance.

Specifications



Shown with 4 way yoke

Product specifications	
Product weight	11kg
Product length	330mm
Product width	220mm
Product height	180mm
Material	Aluminium / ABS
Safe working load	454kg
Colour	White / grey
Warranty	12 months
Power	240V



4 button hand control



Transactive screen



Suits loop slings