

Flexi-Link arm reacher

Efficient, safe and simple to operate. Our patented Flexi-Link connection system is today's most advanced device for connecting the portable lift to the trolley. It's removable wand allows carers to move around freely as there are no handles in place and travel of the hoist is not impeded by lights, fans or doors. Available in 600mm (24") or 900mm (36") length wands.

Accessory is suitable for use with Portable Ceiling Lifts.

Specifications



Product specifications	
Product weight	3kg
Product length	600 / 900mm
Product width	75mm
Product height	30mm
Material	Aluminium /steel
Safe working load	300kg
Colour	Gold / silver
Warranty	12 months



Flexi-link engaged



Flexi-link head