

Trolley table for Luna hoist

Using the Luna ceiling hoist in a portable mode allows for unparalleled flexibility in patient transfer situations, in both nursing homes and home care environments.

The Luna Trolley is a quick, efficient and practical means of moving the portable Luna ceiling hoist from one room to another. With the Luna Trolley, there is no need for the carer to carry the overhead lift from room to room. The smooth rolling Luna Trolley does the job for you.

The Luna Trolley provides the carer with a complete portable solution: a cradle to securely hold the Luna ceiling hoist in place, a convenient on board charger, an extension arm for unhooking the Luna, and a useful pocket for carrying slings or documents.

Specifications



Product specifications	
Product weight	10kg
Product length	370mm
Product width	490mm
Product height	950mm
Material	Powdercoated steel
Colour	Silver / black
Warranty	12 months



Rail dimensions



Straight rail



Luna hoist