

Luna fixed / portable hoist

The Luna overhead tracking hoist is unique in making a track hoist into a flexible technical aid. Luna can be mounted in different positions on the track and may be changed at will to suit varying care needs. The Luna track unit can also be adapted to suit existing tracks so that total replacement is not necessary.



Specifications

Product specifications	
Product weight	7kg
Product length	320mm
Product width	320mm
Product height	130mm
Material	Aluminium / ABS
Safe working load	200 / 275kg
Colour	Silver / black
Warranty	12 months
Power	100 to 240V



Luna Portable



Luna hoist track