

Weigh Scale

Weigh scale for all purposes

Primary application

Home care
Nursing homes
Hospital

The weigh scale is a M weigh (medically approved weight). The weigh scale is designed to weigh disable patients body weight, when they are placed in a sling.

The weigh scale has a big screen and multiple functions including On/Off/Zero/Tare button and Hold button. Weighing capability is up to 300 kg but the LOWEST weight capacity decides the maximum lifting capacity. The scale has an accuracy of 0,1kg.

With the function Hold it is possible to sample 5 to 10 sets of reading and analyse the discrepancy.

The Tare function makes it possible to subtract additional weight like clothing from the total weight.

Key standard functions

On/Off/Zero/Tare button
Hold button
Big display
Weighing capability up to 300 kg
Accuracy of 0,1kg



On/Off/Zero/Tare and Hold buttons




Big display




Easy to use and read

TECHNICAL SPECIFICATION



Weigh scale for Luna with strap




Weigh scale for Luna with bayonet



Weigh scale on Solar 3-point spreader bar



Weigh scale on Solar 2-point spreader bar



Consult the Technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Active Mobility Systems reserves the right to changes without notice.

Product photos may include accessories. Consult with Active Mobility Systems.

