


Manual / Electric Turntable

The Multi-Port Powered Turntable offers the end user ultimate flexibility in situations that require access to a number of transfer locations. The turntable provides 10 take-off points in graduations of 22°. Suitable for hoists up to 300kg SWL. Electric version operated via hoist controls, no need for additional power supply or manual version operated by simple pull cord mechanism.

Specifications



Product specifications	
Product weight capacity	Up to 300kg
Product length	400mm
Product width	400mm
Product weight	5kg
Material	Aluminium / steel
Colour	White
Warranty	1 year against



Rail dimensions



Straight rail



Luna hoist