

Allambie Physiotherapy Treatment Table – 2 Section

This Australian made treatment table exceeds all expectations. With a full size table top it is long enough and wide enough to suit majority of clients. The minimum height of just 50cm makes access for frail, elderly and disabled clients easy and safe. The maximum height of 100cm is high enough for the tallest carers to work ergonomically and comfortably.

4 individual braked wheel system provides a stable and secure platform. One directional locking wheel helps control table when moving.

155mm clearance below frame allows easy and clear access for patient lift hoists.

350kg SWL and 710mm width options and head massage section are available.



Specifications

Product specifications	
Product weight capacity	200 / 350kg
Product length	1930mm
Product width	600 / 710mm
Product height	500 to 1000mm
Product weight	70kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Power	240V
Warranty	5 years parts only



Head massage option



Head section option



Vinyl cover options