

### Avalon Traction Couch – 3 Section

This Australian made traction couch exceeds all expectations. With a full size table top it is long enough and wide enough to suit majority of clients. The minimum height of just 54cm makes access for frail, elderly and disabled clients easy and safe. The maximum height of 104cm is high enough for the tallest carers to work ergonomically and comfortably. The traction unit mount has universal fittings for the major traction units, is height adjustable to improve angle of pull as required and can be swung to either side to improve access.

4 individual braked wheel system provides a stable and secure platform. One directional locking wheel helps control table when moving.

155mm clearance below frame allows easy and clear access for patient lift hoists.

350kg SWL and 710mm width options are available.



#### Specifications

Product specifications	
Product weight capacity	200 / 350kg
Product length	2000mm
Product width	600 / 710mm
Product height	500 to 1000mm
Product weight	85kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Power	240V
Warranty	5 years parts only



Vinyl cover options