

Empire Bariatric Neurological Table

This Australian made bariatric neurological table exceeds all expectations. With a safe working load of 400kg this table will cater for every client. The minimum height of just 50cm makes access for all clients easy and safe. The maximum height of 90cm is high enough for the tallest carers to work ergonomically and comfortably. Electric height foot switch adjustment allows ease of use while caring for client. This table also offers $\pm 15^\circ$ trendelenberg action.

155mm clearance below frame allows easy and clear access for patient lift hoists

Optional hand control or both foot and hand control, 1 or 3 breathe holes with plugs and 100, 110 or 120cm widths is available



Specifications

Product specifications	
Product weight capacity	400kg
Product length	2000mm
Product width	1000, 1100, 1200mm
Product height	500 to 900mm
Product weight	150kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Warranty	5 years parts only



Vinyl cover options