

Enfield Neurological Table – 3 section

This Australian made neurological table exceeds all expectations. With a full size table top it is long enough and wide enough to suit majority of clients. The minimum height of just 45cm makes access for frail, elderly and disabled clients easy and safe. The maximum height of 95cm is high enough for the tallest carers to work ergonomically and comfortably. Electric height foot switch adjustment allows ease of use while caring for client. The electric adjustable head and foot section allows easy positioning of client.

Available in 1000mm, 1100mm or 1200mm widths.


155mm clearance below frame allows easy and clear access for patient lift hoists

Optional hand control or both foot and hand control, 350kg bariatric model, and locking castors are available



Specifications

Product specifications	
Product weight capacity	200 / 350kg
Product length	1930mm
Product width	1000, 1100 , 1200mm
Product height	450 to 950mm
Product weight	85kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Warranty	5 years parts only



Vinyl cover options