

Kensington Adult Change Table

This Australian made change table exceeds all expectations. With a full size table top it is long enough and wide enough to suit majority of clients. The minimum height of just 45cm makes access for frail, elderly and disabled clients easy and safe. The maximum height of 95cm is high enough for the tallest carers to work ergonomically and comfortably.

An 800mm width, and 350kg model is also available

4 individual braked wheel system provides a stable and secure platform. One directional locking wheel helps control table when moving.

155mm clearance below frame allows easy and clear access for patient lift hoists

Additional Cot side and Cot end options available



Specifications

Product specifications	
Product weight capacity	200/ 350kg
Product length	1930mm
Product width	710mm
Product height	450 to 950mm
Product weight	70kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Warranty	5 years parts only



Cot side option



Cot end option



Vinyl cover options