

Leichhardt Junior Change Table

This Australian made change table exceeds all expectations. The minimum height of just 45cm makes access for frail and disabled clients easy and safe. The maximum height of 95cm is high enough for the tallest carers to work ergonomically and comfortably.

4 individual braked wheel system provides a stable and secure platform. One directional locking wheel helps control table when moving. Foot controlled adjustment makes it easy and safe for the carer to adjust without bending.

155mm clearance below frame allows easy and clear access for patient lift hoists

Additional Cot side and Cot end options available



Specifications

Product specifications	
Product weight capacity	200/ 350kg
Product length	1250mm
Product width	710mm
Product height	450 to 950mm
Product weight	70kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Warranty	5 years parts only



Cot side and end options



Swingaway sides



Vinyl cover options