

SpaceSaver

The height adjustable SpaceSaver change table has been designed and manufactured in Australia as an alternative to traditional wall mounted units.

This freestanding unit removes the need for specialised installation and makes shifting and cleaning easy. It compacts to less than 30cm off the wall when not in use which is ideal where space is already compromised.

Patient transfer and changing is simplified with the 450 to 850 mm height range.



Specifications

Product specifications	
Product weight capacity	150kg
Product length	1200, 1500 or 1800mm
Product width	700mm
Product weight	100kg
Lifting Range	450 – 850mm
Control options	Foot pedal, handset or cordless remote
Upholstery	High quality medical grade vinyl over 32mm foam
Frame	Powdercoated steel
Warranty	5 years excluding upholstery



Spacesaver in raised position



Spacesaver with platform in storage position



Optional padded cot cot sides