

Surrey Examination Couch – 2 section

This Australian made examination couch exceeds all expectations. With a full size table top it is long enough and wide enough to suit majority of clients. The minimum height of just 45cm makes access for frail, elderly and disabled clients easy and safe. The maximum height of 95cm is high enough for the tallest carers to work ergonomically and comfortably. Manual or electric height adjustment allows ease of use while caring for client. The adjustable head section allows easy positioning of client. Available in 600mm or 710mm widths and 200 or 350kg safe working load. A range of accessories are available including towel roll holder, cot sides, treatment arm boards, height adjustable armrests and knee supports.



Specifications

Product specifications	
Product weight capacity	200 / 350kg
Product length	1930mm
Product width	600 or 710mm
Product height	450 to 950mm
Product weight	70kg
Material	Powder coated steel / foam
Frame colour	Light grey
Vinyl colour	Std blue – others available
Warranty	5 years parts only



Vinyl cover options