

Comet Ultra

**BUILT TO SUPPORT WITH PERFORMANCE
AND RELIABILITY**



Ideal for those with larger body shapes, the Comet® Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality. The Comet Ultra will get you to your destination quickly and efficiently, in complete comfort and control.



Active Mobility

Comet Ultra



The safety and stability of the Comet Ultra is never a concern; even for users weighing up to 220kg. Developed with the highest quality components and safety features, the Comet Ultra makes a trustworthy companion no matter where you go. Furthermore, its flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference.

The special contoured Ultra Seat has been developed in conjunction with therapists to ensure an optimum seating position and increased comfort. In addition, a pneumatically adjustable lumbar support allows for further comfort when seated if and when required. The reinforced seat base and adjustable armrests offer excellent support and allows the use of the armrest as additional support when getting on and off the scooter.

Features



Improved suspension to offer superior comfort to larger individuals.



















Reinforced seat column offers additional stability and safety when seated.



Reinforced height and width adjustable armrest.

Tech Specs

 Backrest width 660 mm	 Seat depth 535 mm	 Floor to seat height Min. 440 mm Max. 510 mm	 Backrest height 515 mm	 Driving unit width 665 - 685 mm	 Total length 1450 mm	 Empty weight 148 kg	 Max user weight 220 kg
 Battery capacity 2 x 12 V/75 Ah	 Turning diameter 2750 mm	 Max climbable obstacle height 100 mm	 Ground clearance 100 mm	 Max safe slope 10°	 Max driving distance 49 km	 Max driving speed 10 km/h	 Motor capacity 650 W / 1800 W

Active Mobility Systems
Unit 8, 110-120 Silverwater Rd
SILVERWATER NSW 2128

Ph: 02 9649 2111
F: 02 9649 8506
E: sales@activemobility.com.au
www.activemobility.com.au

