

### Crystal Heavy-Duty Shower Trolley

The Crystal power operated heavy-duty shower trolley was specifically designed to accommodate the larger heavier patient. At 2000 mm long, it is 100 mm longer than the standard Crystal and it boasts an adjustable mattress platform width of 980 mm. The safe working load has also been substantially increased to 205kg. This product has a fixed angle mattress for effective patient positioning and a vertical flip function for space saving and easy cleaning access.

#### Specifications

The Crystal heavy-duty shower trolley is equipped with lockable fold-down side rails for optimum safety, dual braking feature and large easy roll castors making it easy to manoeuvre



Product specifications	
Product weight capacity	205kg
Product length	2200mm
Product width	1100mm
Product weight	70kg
Lifting range	590 – 1000mm
Control options	Hand control or foot control
Upholstery	Waterproof removable mattress in soft PVC
Frame	Steel frame, grey powdercoated, anti corrosion treated
Warranty	12 months

