

## Crystal Paediatric Shower Trolley

The 1600mm Crystal Paediatric Shower Trolley combines speed and silence. It features a mattress designed for optimised washing and an adjustable headrest.



### Specifications

Product specifications	
SWL	160kg
Weight	65kg
Length	1600mm
Width	780mm
Height	560 – 1000mm
Tilt Range	0 – 6.5o
Warranty	12 months

### Features

- The patient's comfort is ensured by the texture of the mattress
- The triangular headrest can be adjusted to any position
- Stainless steel side rails ensures the patient's safety