

Disposable Kwik-Fix Slings

Single patient use sling can be used as a general purpose sling for short stay at hospital or residency visits, supplied with a stiffened head support. Easy to apply, all supplied with padded legs and additional tailoring to commode aperture as standard. Used specifically in acute markets to limit cross infection.

Specifications



Product specifications	
Transfer type	General transfer
Lift type	Active & Passive
Support type	Head & Torso
Attachment type	Loop fixing: Loop fixing spreader bar & Clip fixing: Stud fixing spreader bar
User physicality	Limited body control, head control, trunk and hip control & mobility
Sizing	Small – Xlarge
Weigh limit	200kg (31 stone/ 440lbs)
Mobility FIM Measure	1
Available material	Coated paper

Sizing Guide

	A (Shoulder Circumference)	B (Top of Head to Base of Spine)
Small	76 - 86cm	69 - 77cm
Medium	86 - 96cm	77 - 85cm
Large	96 - 108cm	85 - 92cm
XLarge	108 - 130cm	92 - 98cm