

Falls Management



Active Mobility



Fall Prevention Alarms



Slippers




Hip Protectors



DM3

ALARMS



Posey Sitter Elite

- Unit alerts staff to patients attempts to exit bed, chair or toilet seat unattended
- Simple to operate and compatible with Posey chair, bed or toilet
- Personalized voice recording capability
- Contains pre-recorded voice message
- Five alarm tones
- Option available for tone and / or voice alarm
- Adjustable volume
- Battery operated with 4 'AA' batteries
- Adjustable cord (790mm - 1520mm) with locking clip (sold separately)
- Over-moulded impact absorbing case
- SUSPEND allows extra time for administering care without alarm activation
- Sensor purchased separately



Posey Keepsafe Deluxe

- Unit alerts staff to patients attempts to exit bed, chair or toilet seat unattended
- A lightweight compact fall alarm
- Compatible with all Posey bed, chair, toilet, floor and PIR sensors
- Five alarm tones
- Battery operated with 4 'AA' batteries
- Over-moulded impact absorbing case
- SUSPEND allows extra time for administering care without alarm activation
- HOLD allows resident to be away from the alarm without having to turn it off



Posey Keepsafe Essential

- Unit alerts staff to patients attempts to exit bed, chair or toilet seat unattended
- A magnet and pull-cord attaches to patients clothing
- Can be used on bed, chair or toilet
- Adjustable cord (790mm - 1520mm) with locking clip
- Option for tone and/or voice alarm
- Five alarm tones
- Pre-recorded voice message
- Adjustable volume
- Battery operated with 4 'AA' batteries
- SUSPEND allows extra time for administering care without alarm activation

ALARM SENSORS

Chair Pad - Sensor



- For patients at risk of falling
 - For wheelchair or other seating use
 - Compatible with Keepsafe Alarm, and Keepsafe Deluxe
-
- Width: 330mm
 - Depth: 330mm

Posey Alarm - Bed Sensor



- Compatible with Keepsafe Deluxe and Elite alarm units
- Quality, reliable sensor to reduce false alarms
- Sensor placed on top of mattress

Posey Alarm - Floor Sensor Mat



- Available in two sizes
 - Small for usage in room or facility door. Cord measures 3050mm
 - Large to be used at bedside. Cord measures 3050mm
 - Compatible with the Keepsafe Deluxe and Elite alarm units
-
- | | | |
|--------|-------|--------|
| • Size | Small | Large |
| Width | 910mm | 1520mm |
| Length | 610mm | 910mm |

ALARMS

Alarm Monitor - Infrared Wireless



- Monitor Detector
- Can be used with a bed or door
- Can be set to 'chime' at local alert or can send a remote signal to the receiver
- Receiver range up to 30m depending on building construction
- Motion detector & receiver is battery operated via 9 Volt battery (battery sold separately)

Wireless Personal Pager



- Push button transmitter / pendant
- Receiver can be set to 'chime' or 'vibrate'
- Receiver range up to 30m depending on building construction
- Battery operated - Transmitter takes 12 Volt
Receiver takes 2'AA' batteries (sold separately)

Toilet Seat Sensor



- For patients at risk of falls
- Compatible with Keepsafe alarm and Keepsafe deluxe
- Placed under toilet seat
- Activates when patient leaves toilet

FLOOR MATS

Floor Cushion - Economy



- For patients at high risk of falls
- Tested for shock absorbency
- Helps absorb up to 65% of the impact from a fall
- High density, polyurethane foam
- Folds in half with convenient carry handle
- Meets CA#117 section E Flame Retardancy Standards
- Size 1830mm x 610mm x 40mm (W x D x H)
- Folded 610mm x 910mm x 80mm (W x D x H)

Fall Mat - Bevelled, High Vis Floor Cushion



- To be used with patients at high risk for falls
- Helps reduce impact forces by up to 85%
- EVA foam
- Wipe clean vinyl cover
- Non-skid bottom
- Tri-fold design with carry handle
- 25mm thick foam creates low profile edge
- Added hard plastic on the bevelled edge helps floor cushion lay flat & reduce curling
- Bevelled edge has glow strip for high visibility
- Meets CA#117 Section E Flame Retardancy Standards
- Folded 610mm x 970mm x 80mm (W x D x H)

Fall Mat - Bevelled Floor Cushion



- To be used with patients at high risk for falls
- Helps reduce impact forces by up to 85%
- EVA foam
- Wipe clean vinyl cover
- Non-skid bottom
- Tri-fold design with carry handle
- 25mm thick foam creates low profile edge
- Meets CA#117 Section E Flame Retardancy Standards
- Folded 610mm x 970mm x 80mm (W x D x H)

SLIPPERS



Non-skid Slippers

- Non-intrusive colours
- Recommended use for patients that require non-skid footwear
- Colour coded by size
- Size available: Small, Medium, Medium/Large, Large, Extra large
- Soft terry cloth upper
- Latex free rubber sole
- Machine washable
- Supplied in pairs

Specifications

Size	Women shoes size	Men shoes size	Colour
Small	4-5 AUS (6-7 US)	N/A	White
Medium	5-7 AUS (8-9 US)	4-7 AUS (5-7 US)	Blue
Medium/Large	8-10 AUS (9-11 US)	8-10 AUS (8-10 US)	Brown
Large	N/A	10-11 AUS (US)	Beige
Extra Large	N/A	12-13 AUS (US)	Navy



Non-skid Slippers - Wide Opening

- Non-intrusive colours
- Recommended use for patients that require non-skid footwear
- Colour coded by size
- Size available: Small, Medium, Medium/Large, Large, Extra large
- Larger opening to help accommodate patients with oedema or large ankles
- Soft terry cloth upper
- Latex free rubber sole
- Machine washable
- Supplied in pairs

Specifications

Size	Women shoes size	Men shoes size	Colour
Small	4-5 AUS (6-7 US)	N/A	White
Medium	5-7 AUS (8-9 US)	4-7 AUS (5-7 US)	Blue
Medium/Large	8-10 AUS (9-11 US)	8-10 AUS (8-10 US)	Brown
Large	N/A	10-11 AUS (US)	Beige
Extra Large	N/A	12-13 AUS (US)	Navy

SLIPPERS



Non-skid Slippers - Wet/Dry, Shower/Pool

- Recommended use for patients requiring footwear in aquatic rehab pools and floors that have any texture
- Ideal for shower or aquatic pool use
- Loose knit nylon lets water run through
- Nylon mesh upper dries quickly when wet and will not hold water
- Stretchable nylon uppers
- Latex free rubber sole
- Non-skid protection when dry
- Supplied in pairs

Specifications

Size	Women shoes size	Men shoes size
Small	4-5 AUS (6-7 US)	N/A
Medium	5-7 AUS (7-8 US)	4-7 AUS (5-7 US)
Large	8-10 AUS (9-11 US)	8-10 AUS (US)
Extra Large	N/A	10-11 AUS (US)
XX Large	N/A	12-13 AUS (US)



Non-skid Slippers - Falls Management Colours

- Colours are easily incorporated into any fall management system
- Clearly identifies to staff patients at high fall risk
- Recommended use for patients that require non-skid footwear
- Available in 5 colours
- Sizes available: Small, Medium, Medium/Large, Large, Extra Large
- Soft terry cloth upper
- Latex free rubber sole
- Machine washable
- Supplied in pairs

Specifications

Size	Women shoes size	Men shoes size	Colour
Small	3-4 AUS (5-6 US)	N/A	Green, Orange, Purple, Red, Yellow
Medium	4-5 AUS (6-7 US)	4-5 AUS (5-6 US)	Green, Orange, Purple, Red, Yellow
Medium/Large	5-7 AUS (7-8 US)	5-7 AUS (6-7 US)	Green, Orange, Purple, Red, Yellow
Large	8-10 AUS (9-11 US)	7-9 AUS (US)	Green, Orange, Purple, Red, Yellow
Extra Large	N/A	10-12 AUS (US)	Green, Orange, Purple, Red, Yellow

FOOTWEAR & FALLS KIT



Non-slip Socks

- Easily incorporated into any fall management system
- Non-intrusive colours assist staff members to clearly identifying patients at fall risk
- Recommended use for patients requiring non-skid footwear
- Available in 5 colours - Green, Orange, Purple, Red & Yellow
- Soft terry cloth
- Latex free
- Washable
- Skid resistant rubber tread on both sides
- Large size features a mesh band that fits large width ankles and calves up to adult size 12



Fall Management Kit

- Kit used for easy identification of high risk fall patients
- All components for high fall risk patients in one pre-packaged kit
- Kit includes blanket, socks, magnet and bracelet
- Available in 5 colours - Green, Orange, Purple, Red and Yellow



Toilet Seat Sensor

- High impact- absorbing viscoelastic pads help protect hip bones against injury from falls
- Soft, comfortable pads
- 100% latex free
- Discreet and low profile are virtually undetectable under clothing
- Washable
- Same look as regular briefs
- Soft pink cotton blend material with decorative lace
- Reinforced panty line
- Removable pads
- Available in 5 sizes: Small, Medium, Large, Extra Large, XX Large

HIP PROTECTOR

Hip Protector - Community Men



- High impact - absorbing viscoelastic pads help protect hip bones against injury from falls
- Soft, comfortable pads
- 100% latex free
- Discreet and low profile are virtually undetectable under clothing
- Washable
- Same look as regular briefs
- Soft grey cotton blend material
- Removable pads
- Available in 5 sizes - Small, Medium, Large, Extra Large, XX Large

Hip Protector - Ez-On



- High impact - absorbing viscoelastic pads help protect hip bones against injury from falls
- Soft, comfortable pads
- 100% latex free
- Discreet and low profile are virtually undetectable under clothing
- Washable
- User can wear their own undergarments
- Can be worn in the shower
- Removable pads
- Available in 5 sizes - Small, Medium, Large, Extra Large, XX Large

Hip Protector - Standard



- High impact - absorbing viscoelastic pads help protect hip bones against injury from falls
- Soft, comfortable pads
- 100% latex free
- Discreet and low profile are virtually undetectable under clothing
- Washable
- Easily fits over undergarment or can be worn as underwear
- Removable pads
- Available in 3 sizes - Small, Medium, Large

LIFTING CUSHIONS



CAMEL Lifting Cushion

- Suitable for very large or bariatric people (SWL 320kg)
- Lightweight
- Simple to use
- Appropriate for any age
- Portable
- Packs away into a compact and discreet bag (carry bag provided)
- For use indoor or outdoor
- Minimises moving & handling risks
- Enhances independent living
- Battery operated



ELK Lifting Cushion

- Designed to be used with the help of a carer the ELK emergency lifting cushion can be used in the tightest of spaces
- Lifts up to 450 kg
- Portable - weighs just 3.6 kg and folds up for easy carrying
- Applicable for any age
- Simple to use
- Minimises moving and handling risk - protects carers' backs
- Suitable for indoor our outdoor use
- Operated by low air pressure using the Airflo power source

Course DM3

Fall Rehabilitation

Course CM3 is designed for the elderly and / or weakened patient.

By offering various modes of ambulation, carried out in a safe controlled setting, Course DM3 stimulates locomotor function, improves balance and preserves self-confidence.

As part of an extended care plan, it is an objective tool for falls prevention, risk evaluation and 'better balance' workshop.

Course DM3 is an ideal investment for aged care facilities and incommunity accomodation for the disabled. It can also be used in the acute hospital rehabilitation setting and occasional day care/respice.



Specifications

- Stimulating locomotor functionality
- Improving balance
- Preserving self-confidence
- Maintaining independence

Evaluation tool:

The evaluation grid provided with the DM3 makes it possible to analyse the progress of each user with reference to the level of assistance required and the degree of success with which each exercise is completed. It also facilitates setting up groups with similar skill levels, thereby promoting social interaction and encouraging friendly competition.

12 exercise based around:

- Physiological walk
- Static balance
- Dynamic balance

Name	CARPET #1. (180 points)				CARPET #2. (54 points)					CARPET #3. (66 points)						Total
	Forward Walk	Right Side Walk	Left Side Walk	Backward Walk	Collecting	Climbing	Descending	Balance on Right	Balance on Left	Transition	Right cant	Left cant	Bars	Carpet	Mushroom	
	Score /45	Score /45	Score /45	Score /45	Score /18	Score /9	Score /9	Score /9	Score /9	Score /18	Score /6	Score /6	Score /9	Score /9	Score /18	
																/300
																/300
																/300
																/300
																/300
																/300

Memo: Multiplying coefficient of 3 if achieved without assistance | Multiplying coefficient of 2 if achieved with technical assistance (canes, or walking frame) | Multiplying coefficient of 1 if achieved with manual assistance of accompanying carer.



Active Mobility

Active Mobility Systems

Unit 8, 110-120 Silverwater Rd., SILVERWATER NSW 2128

P: 02 9649 2111 F: 02 9649 8506

W: www.activemobility.com.au E: sales@activemobility.com.au

