

Iris



Variable Height Profiling Shower Table

Iris - is a new and innovative product from Reval. Iris has been designed to facilitate shower bathing and hygiene routines, for patients and service users with a wide range of postural and cognitive challenges. Its three section profiling mattress platform can be adjusted to provide specific postural support in the four key regions of the anatomy, head, back, seat, and legs providing maximum safety, support, and comfort. Iris is also kind on the care provider thanks to its variable seat and platform height adjustment. Care workers can maximise their comfort and safety using this function which enhances their wellbeing. Iris is also a mobile patient transfer system and can be used to transfer patients and service users in the horizontal or seated position. Reval manufactures a wall mounted shower and sluice systems. When incorporated into the room design, these appliances help to maximise the working space, creating a more practical working environment.



Iris




Shower system



Shower and sluice system

Variable Height Profiling Shower Table



SPECIFICATIONS:

Chassis

- Aluminium & paint finished rolling chassis.
- Aluminium actuator column.
- ABS plastic chassis covers.
- x 4 Guiding and braking castors.
- Steering handle.

Platform

- x3 Articulating patient support sections in polypropylene.
- Stainless steel articulating hinges and fixations.
- x6 Stainless steel security side rails with locking functions.
- Stainless steel foot end rail.
- Tilting stainless steel head rail.

Mattress

- Soft PVC waterproof mattress designed for heavy usage.
- Abdominal safety belt.
- Drainage outlet and waste water hose at foot end Ø40mm.
- x2 Triangular headrests, in soft durable waterproof PVC material.

Electric components

- 24Vdc IPX6 rated variable height platform actuator
- 24Vdc IPX6 rated backrest actuator
- 24Vdc IPX6 rated legrest actuator
- 24Vdc IPX6 rated retracting chassis actuator
- IPX6 Control module
- IPX6 8 function hand set
- 5Ah / 24V Battery
- 24V battery and charger module.

Performance:

- SWL – Max patient weight 160Kg
- Weight of appliance – 74Kg
- Usage capacity - max 35 cycles @ 80Kg

Conformity:

- 93/42/CEE (Medical Equipment Directive)
- NF EN 12182 (Tech. support systems for less able persons)
- NF EN 1441 (Risk analysis for medical equipment)

Warranty:

- 2 years.



Bathing and Hygiene



- Multi-position posture management
- Variable height mattress platform
- Comfortable profiling mattress
- Waterproof mattress with high sides
- Lockable security side rails
- Comfort pillow and body belt
- Retracting stabilising chassiss

Variable Height Profiling Shower Table



Stay with us...

For more information, please call or email AMS team

Active Mobility Systems

Freecall 1800 543 837

Email sales@activemobility.com.au

www.activemobility.com.au



Unit 8, 110-120 Silverwater Rd, SILVERWATER NSW 2128

info@activemobility.com.au, www.activemobility.com.au