

Jacaranda Bed



Active Mobility

Features

- Electric adjustable bed ideal for homecare
- Electric high/low functions, including backrest, knee-bend and vertical lift
- Additional low-low feature reduces the risk of falls
- Electric Trendelenburg and reverse Trendelburg positions
- Pull-apart function for easy transport and installation
- Supplied with European Beech head and footboard panels
- Hand control is simple and easy to use
- 4x 125mm locking castors
- Battery back-up
- Beige powder coated steel frame

Specifications

Overall Width: 900mm
Overall Length: 2080mm
Sleeping Surface: 1980x870mm
Height: 200 - 765mm
SWL: 180kg

Accessories Available

Horizontal Side Rails (Set)
Self Help Pole
Clamp on Rail

