



# TILT IN SPACE REHAB SHOWER COMMUNE

Improving your  
life everyday™

VERSION 1.0 09/2020



Active Mobility

# TILT IN SPACE

## REHAB SHOWER COMMUNE



Tilt in Space Commode model shown with optional extras - seat, laterals, knee supports, armrest pads, pommel/splash guard. Options sold separately.

Lateral thoracic support adjustable height mounting plates standard on all models



## Features

- 445mm and 500mm seat widths
- Lateral thoracic support adjustable height mounting plates standard on all models. 240mm range (20mm increments)
- Height adjustability of the frame
- Knee guide mounting plates standard on all models. 40mm horizontal mounting range (20mm increments)
- Height adjustable flip up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Standard left and right side opening on seat frame
- Durable and long lasting with a stainless steel frame
- Swing away footrests with padded calf supports
- 4 x 125mm locking swivel castors



Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Frame Height Range	Max User Weight
Adult Tilt in Space Rehab Commode <b>KA11507-REHAB</b>	445mm	435mm - 465mm	855mm	580mm	1070mm	540mm - 465mm	200kg
Adult Tilt in Space Rehab Commode <b>KA11512-REHAB</b>	500mm	435mm - 465mm	855mm	635mm	1070mm	540mm - 465mm	200kg

Conforms to AS/NZS 3973:2009 for safety and peace of mind



Active Mobility

**Active Mobility Systems**

Unit 8/110-120 Silverwater Road, Silverwater NSW 2128

**Ph:** 02 9649 2111

**Fax:** 02 9649 8506

**Email:** [sales@activemobility.com.au](mailto:sales@activemobility.com.au)

**Web:** [www.activemobility.com.au](http://www.activemobility.com.au)

**Improving your  
life everyday™**