

Lambda Shower Trolley

The Lambda Shower Trolley is designed to help bedridden clients go to the bathroom while at the same time provide a good working environment for carers when heavy lifting must be eliminated. It is extremely easy to use since all the functions can be found on the hand control. A separate motor allows the showering surface to tilt up to 10o thus ensuring client comfort and adequate water drainage



Specifications

Product specifications	
SWL	175kg
Weight	115kg
Length	2060mm
Width	780mm
Wheels	150mm
Material	Powder coated steel
Warranty	12 months

Features

- Hi-Lo battery powered lifting
- Anti-Trendelenburg as standard
- Adjustable and lockable side and end rails
- All functions controlled by hand control
- Large and easy to manoeuvre wheels