

4 Point Yoke

The Allegro 4 point yoke is a high quality attachment for a floor hoist which allows the carer greater control over moving and lifting the patient. The 4 point yoke allows more room and increased comfort for the patient when hoisting. The 4 point yoke is suitable for use with all loop slings.



Specifications

Product specifications	
Product weight capacity	200kg
Product length	500mm
Product width	350mm
Product height	170mm
Product weight	3kg
Material	Powder coated steel / foam
Colour	Light grey / green
Cover	Medical grade vinyl
Warranty	1 Year excluding cover



Fits to floor hoists



Removable vinyl cover



Suits loop slings