

Pivot Frame

The Allegro pivot frame is a high quality attachment to a floor hoist which allows the carer to easily move the patient to a sitting or lying position with minimum effort. This pivot frame uses slings with the keyhole clips for greater safety and ease of connecting the sling to the frame

Specifications



Product specifications	
Product weight capacity	200kg
Product length	575mm
Product width	700mm
Product height	500mm
Product weight	3kg
Material	Powder coated steel
Colour	Light grey
Cover	Self skinning polyurethane
Warranty	1 year excluding upholstery



Fits to floor hoists



Keyhole buckle



Suits keyhole buckle slings