

SONATA 150

Compact Mobile Patient Lifting Hoist



Innovation in
Patient Handling

Allegro
CONCEPTS PTY LTD

SONATA 150

Sonata 150: A lightweight compact mobile patient lifting hoist with a safe working load of 150kg and an innovative and ergonomic design.

The Sonata 150 is an easy to use everyday lifter providing superior patient lift space through the use of SPACE LINK technology.



SPACE LINK™

- breakthrough for optimized patient space
- enhances floor lift capacity

Sonata 150 Specifications

Safe Working Load	150kg
Base Length	1005mm
External Base Width Closed	605mm
External Base Width Open	1070mm
Internal Base Width Closed	510mm
Internal Base Width Open	975mm
Base Width Adjustment	3 Position Manual
Under Bed Clearance	130mm
Minimum Height Yoke	420mm
Maximum Height Yoke	1735mm
Average Lift Range*	870mm
Lifter Weight	41kg
Minimum Storage Height	1260mm
Under Base Height	55mm
Maximum Hoist Reach	875mm
Minimum Hoist Reach	555mm
Hoist Reach at Maximum Height	555mm
Mass of Heaviest Part	22kg (Mast assembly)

* With general purpose large sling.

The Allegro Concepts range of patient lifting equipment including lifters, slings and accessories has been designed and tested to comply with: AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.



Standard Yoke




4 Point Yoke



Tri Hook Yoke



Hilift Yoke



Pivot Frame



Weigh Scale

Allegro Concepts Pty Ltd

Unit 9, 110-120 Silverwater Road Silverwater NSW 2128

Phone: (02) 9749 7812 Fax: (02) 9749 2144

Email: sales@allegroconcepts.com.au

Website: www.allegroconcepts.com.au

Allegro
CONCEPTS PTY LTD
www.allegroconcepts.com.au