



# HEADREST SUPPORTS

(Please store this in your Medifab Folder)

Product Code 1270

1270



- ✔ Lightweight, pressure relieving headrest pads offer comfort and support.
- ✔ Headrest covers are removable for cleaning.
- ✔ Aluminium headrest inners provide strength.

## Multi-Adjustable Headrest

Unique design that doesn't have a protruding mechanism at the rear.



Quick release lever allows the headrest to be removed in an instant.



- ✔ Quick-release
- ✔ Easy to install
- ✔ Strong components
- ✔ Huge adjustment range
- ✔ Fine adjustments for optimal positioning are 'locked-off' securely.

## Basic Headrest System

Basic headrest mechanism with height and depth adjustment only.



## Headrest Pad Options

### Square Pad

Basic, compact headrest for clients with good head control.



### Contoured Pad

Offers close lateral & occipital support for clients with reduced head control.



### Multi-adjustable Lateral Support Pad

Lateral wings of headrest can be adjusted and locked at the desired position.



