



Molift Smart 150

Light, foldable and transportable without tools

Compact smartness

Molift Smart 150 is designed to be easily wheeled and transported to where it is required. It is designed for home care, for travelling and easy storage. The Molift Smart is foldable in one piece without using any tools.

Due to its low weight, the small footprint and the unique construction, Molift Smart 150 is easy to manoeuvre even in narrow spaces. The hoisting range of 270–1680 mm provides an excellent maximum hoisting height and at the same time gets very low, facilitating hoisting from the floor. The 4-point sling bar ensures a comfortable and spacious hoisting position for the user.

Molift Smart 150 has been awarded for Design Excellence by the Norwegian Design Council.



Product name	Item no.	Dimensions		
Molift Smart 150	1020	Length	1200 mm	
		Leg height	110 mm	
		Dimension folded	L1160 x W475 x H360 mm	
Material		Hoisting height	270 - 1665 mm	
Aluminium, steel		Hoisting interval	1410 mm	
Weight		Technical		Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85°C for 3 minutes.
Total Weight	26 kg	Service Software	Included	
Weight when separated	Part 1: 13 kg Part 2: 13 kg	Battery	NiMH 14.4 V - 2.6 Ah	
		Battery charge time	approx. 3 hours	
		Hoisting speed	80 mm/sec	
		Rated performance	50 hoists/75 kg/500 mm	

Accessories



Molift Soft Travel Bag

The Molift Soft Travel Bag protects and enables easy transporting of the hoist, battery, charger and sling during travel.

Item no: 3049177



Travel Suitcase

This hard shell suitcase provides additional protection when transporting by air. It is TSA approved for air travel and ideal for the client with an active life.

Item no: 0990200

Features



4-point sling bar

The standard 4-point sling bar ensures a comfortable and spacious hoisting position for the user.



Foldable without tools

Molift Smart 150 is designed to be easy and convenient at home and to transport when travelling. The possibility to easily fold the hoist without tools is a big advantage.



Active Mobility Systems
Unit 8, 110 - 120 Silverwater Road
SILVERWATER NSW 2128

Email: sales@activemobility.com.au
Phone: 02 9649 2111



Active Mobility