

Monarch Ceiling Hoists by Active Mobility Systems



Active Mobility



Active Mobility




Monarch Fixed Ceiling Hoist

The Monarch Fixed Ceiling Hoist is a complementary product to other products for individuals with mobility challenges. It is designed for use by caregivers and health care professionals to lift and transfer a person in a homecare, assisted-living or hospital environment.

The Monarch Fixed Ceiling Hoist is one of the only ceiling hoists on the market using reliable lithium ion (Li-ion) battery technology. Most competitor products use Lead-acid (heavy, toxic, low cycle life) or NiMh.

Unique Features

- Unique Spreader Bar with Secure Sling Mechanism
- Fast – lifts up to 30% more quickly than other models
- Low lift clearance – slim profile means more room to maneuver patient
- Li-ion battery – recharges twice as fast as average competitor products.
- Ultra Lightweight – Weighing 8kg 35% lighter than competitors
- 200kg safe working load
- Ideal for Nursing Homes, Hospitals or Domestic Use



Fast Charging Li-ion battery



200kg



Low lift clearance



Secure Sling Mechanism



Easy-to-use controls




Magnetic Handpiece

Product Specifications Monarch Fixed Hoist

Specifications

Model	Monarch Fixed Ceiling Hoist
Function	Patient Lifting, Room-to-Room Transfer
Total Weight	7.5 kg (Manual) 8.5kg (Motorised)
Weight Capacity	200 kg Standard
Lifts per Charge	30 to 40 cycles
Speed with 200kg	4.5 cm/s
Soft Start / Stop	Yes
Emergency Stop / Lower	Yes
Track Options	H system, Turntable, Gate, Curves, Room-to-Room, Freestanding Gantry : 2m, 3m, 4.5m, 6m, 8m
Track Compability	Compatible with all Monarch track and other leading track manufacturers
Application	Nursing Home, Hospital, Domestic, School, Leisure Centre, Physiotherapy

Dimensions



Active Mobility



Active Mobility



Monarch Portable Ceiling Hoist

The Monarch Portable Ceiling Hoist is one of the lightest on the market weighing only 4.8kg, whilst lifting the industry standard of 200kg. Ideal for occasional or multiple room use, the product is robust and charges with an AC Charger. The unique spreader bar includes our industry leading Secure Sling Mechanism designed to ensure loop attachments cannot be unintentionally removed.

Perfect for rental properties, multi-room, temporary or short-term use, the Monarch Portable Ceiling Hoist can be used as a flexible alternative to fixed overhead hoists.

Unique Features

- Ultra-Lightweight, Weighing only 4.8kg
- Increased Battery Capacity
- Carry handle
- Express Recharge
- Unique Spreader Bar with Secure Sling Mechanism
- 200kg safe working load
- Ideal for transporting with the user
- Ideal for rental properties or hotels for nonpermanent use




Easy-to-use Carabiner clip



Battery is 4x more powerful



50% more lift cycles



Available with telescopic reacher



Easy Transport Handle




Unique Secure Sling Mechanism

Product Specifications Monarch Portable Hoist

Specifications

Model	Monarch Portable Ceiling Hoist
Function	Patient Lifting, Room-to-Room Transfer, Temporary Use, Rental Properties
Total Weight	4.8 kg
Weight Capacity	200 kg Standard
Lifts per Charge	30 to 40 cycles
Speed with 200kg	4.5 cm/s
Travel	2300 mm
Soft Start / Stop	Yes
Emergency Stop / Lower	Yes
Options	Standard Portable Trolley (Monarch track compatible), Bespoke Portable Trolley
Track Options	(compatible with leading track manufacturers), Reacher.
Track Compability	H system, Turntable, Gate, Curves, Room-to-Room, Freestanding Gantry : 2m, 3m, 4.5m, 6m, 8m
Application	Compatible with all Monarch track and other leading track manufacturers Nursing Home, Hospital, Domestic, School, Leisure Centre, Physiotherapy

Dimensions



Active Mobility



Active Mobility



Monarch Fixed Ceiling Hoist Track

The Monarch Fixed Ceiling Hoist Track offers a wide range of Ceiling Track solutions and accessories to allow installation in virtually any setting.

Options include room covering systems (also known as XY-Systems or H-Systems) that can connect up to a neighbouring room with Ceiling Hoist gate system, tight curves for small bathrooms and a flexible turntable with 16 exit possibilities.

Unique Features

- Off-white Colour to Blend into Ceiling
- Multi-Room Gate System
- Unique No-Show Internal Bracket
- Full room covering XY (H) Systems
- Turntable Systems
- Modular Tracking Units
- Bespoke Pieces Available
- Recessed Tracking




Wall Mounted Track



90° Curve Track



Magnetic Gate System



Flush X/Y (H) System



Turntable System



Flush / Recessed Tracking

Product Specifications Monarch Fixed Track

Specifications

Model Wall Mounted or Freestanding Track

Function Patient Lifting, Room-to-Room Transfer, Ideal for Bathrooms, Bedrooms & Home Use

Dimensions 55mm x 90mm, 140mm or 180mm

Weight Capacity 200kg, 250kg or 400kg

Application Hospital, Nursing Home, School, Domestic, Swimming Pool

Fixing Options Wall Posts, Wall Mounted, Ceiling Mounted, Suspended Ceiling, Sloping Ceiling, Freestanding Gantry

Configuration Options




Straight Track




45° Track




90° Track



Wall Bracket



Gate System







Turntable



Wall Post

Dimensions



Active Mobility

CEILING HOIST SYSTEM

Parts and Descriptions

Fixings

- 17105689 Fixing into Solid Concrete with HKD Anchor
- 17105691 Chemical Resin Fixing with mesh sieve & insert for hollow core concrete
- 17105688 L Bracket kit and Channel with Fixing kit for TGI Joists (1 per fixing)
- 17105850 Coach Bolt with thread- 10X140mm Stud Screw

- ### Wall Post
- 17105681 Wall Post
 - 17105686 Wall Post Diagonal
 - 17106883 Wall Post Skirting

- ### Ceiling Plate
- 171021601

- ### Wall Bracket
- 17107001

- ### Track
- | Rail | 2m x 90mm |
|----------|-----------|
| 17105150 | Rail |
| 17105160 | Rail |
| 17105670 | Rail |
| 17105180 | Rail |
| 17105190 | Rail |

- 17105240 End Cap
- 17400790.1 Rail End Stop

- ### TBC
- Charger

- ### Charging Sensor
- 17400200

- ### Infra-red Kit
- 17000132

- ### Inside Trolley Set
- 17000180

- ### Quick Lock Set
- 174001201

- ### Fixed Hoist
- | Monarch Fixed | 200kg 2FB |
|-----------------|-------------------------|
| M92MONFIX2FB | |
| M92MONFIX4FB | Monarch Fixed 200kg 4FB |
| M92MONFIX2FB272 | Monarch Fixed 272kg 2FB |
| M92MONFIX4FB272 | Monarch Fixed 272kg 4FB |

- ### Wegh Scales - (300kg)
- OPE-2500-304

- ### Spreader Bars
- 32000630 X-Shaped Spreader Bar
 - 32000020 Manual Cradle
 - 3210550 Spreader Quick Release

- ### Unistrut
- 17105684 Unistrut (Above Ceiling, Unpainted)
 - 17105683 Unistrut (Below Ceiling, Painted)

- ### Portable Hoist
- 17000133 Monarch Portable Hoist - 130kg
 - 17000134 Monarch Portable Hoist - 200kg
 - 170001601 Portable Trolley
 - 700-08310 Reacher
 - 700-08310 Removable Hoist Trolley
 - 17105210 Regent Gantry Track

- ### Lateral Brace
- 17101770

- ### Visible Brace
- 700.11380

- ### XY Trolley System Set
- 17400830.1

- ### XY Trolley Charger Interface
- 17401820

- ### Gate System
- 17000140


- ### Joint Clamp
- 17102150.1

- ### Spring Pin
- 17102040

- ### Powered Turntable
- 17000150

- ### 90° Curve
- 17105250

- ### 45° Curve
- 17105260



MONARCH GANTRY SYSTEM

Parts and Descriptions






Track

17105150	Rail	2m x 90mm
17105160	Rail	3m x 90mm
17105180	Rail	4.5m x 140mm
17105190	Rail	6m x 180mm





Active Mobility



Part Number

Description

17105150
17105160
17105670

2 Meter (90mm) rail
3 Meter (90mm) rail
4 Meter (90mm) rail

17105679

4.5 Meter (140mm) rail

17105190

6 Meter (180mm) rail

17105260

45° Curve (90mm)


17105250

90° Curve (90mm)

17101770

Lateral Brace

Product Specifications Ceiling Track Elements




Part Number

17000150


Description

Powered Turntable




17101700

Wall Bracket



17101631

Rod Cover




1740800

Charge Sensor



17400790

Rail Joint



17102150

End Stop




17105240

End Cap



17102180

Cover Shim



200-11140




Ceiling Plate



17400120

Quick Lock Set

Ceiling Track Sets

	Part Number	Description
	17000140	Gate System Set
	17000180	Flush XY Trolley System Set
	17400830	XY Trolley System Set

Ceiling Track Wall Posts




Part Number

17106883
17105681
17105686

Description

Wall Post Skirting
Wall Post
Wall Post Diagonal

Dimensions



Standard Height: 2500mm
Other Heights: Made to bespoke size
Depth: 40mm
Width: 60mm



Active Mobility



Active Mobility



Monarch Art Ceiling Hoist Decal

Bringing the fun inside with high-impact customisable designs, Monarch Art is a perfect way to help create a sensory environment, personalise a bedroom or bathroom, making patient transfer equipment fun.

Monarch Track Art has a wide scope for application, perfect for a child's bedroom or the bedroom of a patient suffering from dementia. With a range of standard designs, where Monarch Track Art really stands out, is its ability to be fully personalised.

Monarch Track can be retrospectively fit to any Monarch tracking systems.

Unique Features


- A unique way to brighten your room
- Available in standard design or bespoke design
- Retrospective application available
- Easy application



Create A Sensory Environment



Bespoke Designs



Choice of pattern or colour



Brighten Your Room



Seamless Application To Track




Retrospective Application Available

Monarch Track Art

Product Specifications

Specifications

Model	Monarch Track Art
Function	Bespoke Designed Decal for Monarch Tracking System
Colour Options	Six Standard Designs: Space Scene, Unicorns, Birds, Clouds, Balloons or Footballs. Bespoke designs available upon request.
Track Compability	Compatible with all Monarch track. Can be retrospectively fit.
Application	Schools, Bedrooms, Nursing Home, Hospital, Domestic, Leisure Centre.





Active Mobility

Active Mobility Systems
Unit 8, 110-120 Silverwater Rd
SILVERWATER NSW 2128

T: 02 9649 2111

F: 02 9649 8506

E: sales@activemobility.com.au

W: www.activemobility.com.au



Proud supporter of

