


## Montar Fold-Up Change Table

The Montar Fold-Up Change Table is a wall mounted examination table that is perfect for small spaces. When folded up it needs very little space leaving the examination room free for other uses. It has an adjustable height range of 600mm and is the perfect examination table for patients with disabilities.



### Specifications

| Product specifications |                              |
|------------------------|------------------------------|
| Product weight (SWL)   | 175kg                        |
| Width                  | 350mm (Open) 1030mm (Closed) |
| Length                 | 1880mm                       |
| Sleeping Surface       | 1880 x 700mm                 |
| Warranty               | 12 months                    |

### Features

- Folds up for space saving
- Drop down safety rail
- Remote infra-red handset