

Nomad portable ceiling hoist

The Nomad is a flexible and user-friendly portable ceiling hoist. It is intuitive and easy to operate. Characterised by an innovative design and user friendliness, the Nomad has a unique and attractive design. The curved shape is designed to afford the user's head ample space while the integrated 4-point yoke gives ultimate safety, comfort and dignity during the lift and transfer. The user will feel safe and stable without being closed in. The combination with a compatible sling ensures the user will achieve a correct and comfortable transfer position.

The battery simply charges by placing hand piece into the dock and automatically shuts off when battery is fully charged.



Specifications

Product specifications	
Product weight	7kg
Product length	500mm
Product width	200mm
Product height	329mm
Material	Aluminium / ABS
Safe working load	230 / 255kg
Colour	White / grey / red
Warranty	12 months
Power	240V (for charger)



Trolley to suit



Charging dock



Suits loop slings