

Pride® Lift Chairs

END-TO-END



Petite to Small Medium to Large Medium to Large Medium to Large Medium to Large Medium Medium to Large Small to Medium Small to Medium Medium to Large Medium to Large Medium Medium Heavy Duty

	LC-440 Power Stand Assist ⁴	C1	C5	C6	C-101	LC-101 Euro Leather	LC-107	D30	560	660	805	T-3	SR-358M	SR-525M	LC-358XL
	Power seat lift	3-Position Lift Chair with fully removable upholstery components 163 cm & below (5'4" & below)	3-Position Chaise Lounger with fully removable upholstery components 163 cm - 183 cm (5'4"- 6'0")	Infinite-Position Chaise Lounger with fully removable upholstery components 163 cm - 183 cm (5'4"- 6'0")	3-Position Chaise Lounger 163 cm - 183 cm (5'4" - 6'0")	3-Position Chaise Lounger with European style rounded arm design 163 cm - 183 cm (5'4"- 6'0")	Infinite-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 163 cm - 183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm (5'4") & below	Infinite-Position Chaise Lounger 163 cm (5'4") & below	2-Position Chaise Lounger Wall Hugger that requires only 10 cm (4") of space between the wall and chair 163 cm - 183 cm (5'4"- 6'0")	Infinite-Position Chaise Lounger featuring Trendelenburg position (elevates legs above head for greater relaxation) and zero- gravity position 163 cm - 183 cm (5'4"- 6'0")	3-Position Chaise Lounger 5'4"- 6'0"	Infinite-Position "Zero-gravity" Trendelenburg Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'6"- 6'1"
Weight Capacity	125 kg (275 lbs.)	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	227 kg (500 lbs.)
Back Width	53.5 cm (21")	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")	74 cm (29")	71 cm (28")	71 cm (28")	76 cm (30")	69 cm (27")	71 cm (28")	69 cm (27")	91 cm (36")
Seat-to-Floor	49 cm (19")	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	51 cm (20")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")
Seat Depth	46 cm (18")	51 cm (20")	46 cm - 51 cm (18-20")	53 cm (21")	51 cm (20")	51 cm (20")	53 cm (21")	56 cm (22")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	56 cm (22")
Seat Width	54.5 cm (21.5")	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20")	51 cm (20")	56 cm (22")	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	66 cm (26")
Top of Back to Seat	69 cm (27")	66 cm (26")	66 cm - 71 cm (26-28")	66 cm - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	70 cm (27.5")	66 cm (22")	63.5 cm (25")	63.5 cm (25")	69 cm (27")	68.5 cm (27")	62 cm (24.5")	67 cm (26.5")	74 cm (29")
Back Style	Blown Fibre	Blown Fibre Accented Pillow-back	Blown Fibre Sewn Pillow-back	Blown Fibre Sewn Pillow-back	Blown Fibre Button-back	Blown Fibre Split-back	Blown Fibre Sewn Split-Back	Blown Fibre Split-back	Blown Fibre Pillow Back	Blown Fibre Pillow Back	Blown Fibre Sewn Split-Back	Blown Fibre Mongolian-back	Blown Fibre Button-back	Blown Fibre Split-back	Blown Fibre Button-back
Fully Padded Chaise	No	Yes ^{1,2}	Yes ^{1,2}	Yes ^{1,2}	Yes ^{1,2}	Yes ^{1,2}	Yes	Yes ¹	Yes ¹	Yes ¹	Yes	Yes	Yes	Yes	Yes
Integrated Battery Backup ³	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Recline Positions	N/A	3	3	Infinite	3	3	Infinite	3	3	Infinite	2	Infinite	3	Infinite	3
Standard Arm Covers	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Standard Head Cover	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Colour Options	Cocoa, Fern, Sandal	Arctic Blue, Chocolate, Oatmeal	Arctic Blue, Chocolate, Oatmeal	Arctic Blue, Chocolate, Oatmeal	Cocoa, Fern, Sandal	Leather: Beige, Brown, Black, Red	Pecan, Burgundy, Ocean	Ocean, Burgundy, Spruce, Latte	Ocean, Burgundy, Spruce	Latte	Cocoa, Fern, Sandal	Black (Leather)	Ember, Timber	Ember, Timber	Burgundy (Vinyl), Latte

2-position is a TV recline.

3-position is a full recline.

Infinite-position is a full recline with independent leg rest function.

Trendelenburg position elevates feet higher than the head.

FABRIC SELECTIONS



SUEDE




LEVANTE VINYL



LEATHER



SERTA PERFECT FABRICS



Note: Measurements may vary +/- 2.5 cm (1").
Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Wood Frame: Engineered Furniture Laminate/Hardwood Frame^{*}
^{*} C5 & C6 frames are made of steel, not Laminate/Hardwood.

1. Removable fully padded chaise.
2. All upholstery components on the C1 are fully removable.
3. Batteries are not included.
4. This chair does not recline.

1. Not all fabrics are available for every model. Please see Specifications for details.