




Active Mobility





www.activemobility.com.au

02 9649 2111



## ECO TRAVELLER Shower Commode



The Eco Traveller is the perfect shower commode for travelling with a foldable frame and 120KG SWL. Still has all the features of other Shower Buddy products including a commode system. A lack of space makes bathroom mobility even harder, especially when travelling. That's why the Eco Traveller's lightweight, but strong aluminium frame can easily be folded for storage and transportation. This, along with its 120kg SWL, makes it ideal for trips away from home.

### Features

- Rolls over most standard toilets
- Fixed height arms
- Removable and swing-away, height-adjustable footrests
- No rust aluminium frame
- Stainless steel parts and hardware
- 4x 100mm locking castor wheels
- Commode system
- Tool-free assembly and easy fold down

### SPECIFICATIONS

Chair width over arms	550mm
Chair width over castor wheels	550mm
Chair width inside arms	450mm
Chair length	840mm including footrest
Backrest width	448mm
Backrest height	400mm (cushion to top)
Max toilet height clearance	470mm
Seat cushion	420x420mm
Seat height from floor	530-580mm
Seat cushion opening	175x335mm

ARTG: 325669 | AS/NZ3973:2009

