



Active Mobility




www.activemobility.com.au

02 9649 2111



ERGO Walking Sling



Suitable for gait training with the **Luna ceiling hoist** the ERGO Walking Sling has a soft padding under the arm section and the overall design is to provide users excellent support and comfort.

Features:

- Used for gait training and general rehabilitation purposes for active users
- Designed to support upper body when raising the user during gait training
- Lift to and from sitting or lying positions
- Removable padding on the shoulder straps and on the leg straps
- Elastic back, ensuring a good fit for the user
- Coloured strap on upper back with size and midline mark, which ensure correct placement of the sling

SPECIFICATIONS

| | |
|-------------------------|---|
| Weight Capacity | 180kg |
| Sizing | |
| Belly Dimensions | XXS: 62-75cm XS: 76-84cm S: 85-92cm M: 93-105cm L: 106-114cm XL: 116-130cm |
| Material | Polyester |
| Spreader Bar | 2, 3 and 4 point |

