



Active Mobility





www.activemobility.com.au

02 9649 2111



KINGSDENE ADULT Change Table



The Kingsdene is designed to address a wide range of manual handling issues to benefit both patient and carers. The height adjustment allows for comfortable handling of the patient to the correct positional height during procedures with enhanced safety and care.

The Kingsdene Change Table includes:

Cot Side - Fitted upholstered cot sides, where both sides are collapsible.

Motor - Utilises a high-grade electric actuator that's driven to enable full adjusted product hi-lo

Twin Lift - Strong twin-lift mechanism ensures rigidity and stability during treatment

Castor - Four 5" castors, all with brakes

SPECIFICATIONS

Construction	High quality steel sections, and sheet construction utilizing the latest manufacturing techniques throughout
Upholstery	Manufactured in line with using high quality foam, lining with vinyl/imitation leather without edge stitches, prohibits/anti bacterial, and also fire retardant
Height Adjustment	Using high-grade electric actuator motor Hi-Lo upholstery platform, elevation 525-975mm
Cot Side	Upholstered cot sides where both sides are collapsible
Coating	Epoxy Powder coating finish
Max. Load	200kg
Weight	Approx 63kg
Dimensions	(L)1970 x (W)710 x (H)525-975mm

ARTG 312129

ISO 9001:2015 ISO13485:2016 ISO14001:2015 ISO 45001:2018 CE

