



Active Mobility





www.activemobility.com.au

02 9649 2111



MONTAR - Fold-Up Change Table



Montar Fold-Up Change Table is a wall mounted examination table perfect for small spaces. When folded up, it needs very little space and leaves your examination room free for other uses.

It has an adjustable height range of 600mm, and the upholstery is industrial grade vinyl, so spills and splashes are not a problem. This is the perfect examination table for patients with disabilities.

The Montar change table also offers remote carer to patient assistance by means of a cordless handset.

Features:

- Folds up for space saving
- Drop down safety rail
- Remote Infra red handset
- Suitable for many applications, such as: *First Aid, Disabled Assistance, Sports Injuries, Medical Assesments, Medical Treatment, Examination and Training*

SPECIFICATIONS		
	STANDARD SIZE	JUNIOR SIZE
SKU	3200000	3200002
SWL	200kg	200kg
Open Overall Dims	1880mm x 1030mm	1400x1030mm
Closed Overall Dims	1880mm x 350mm	1400x350mm
Sleeping Surface	1880 x 700mm	1400x700mm
Height Range	950-350mm	950-350mm
Top Length	1880mm	1400mm

The Montar is manufactured to exceed Australian Standards: **IEC60601-2-52**

