




Active Mobility




www.activemobility.com.au

02 9649 2111



PRESSALIT - Nursing Bench



The Pressalit Nursing Bench 3000 series provides for maximum safety and comfort. Besides an award winning design, the series includes comfortable and durable surfaces, including polyurethane foam and aluminium. The bench is electrically height to suit individual user requirements. Gas assisted tilt allows carers to easily fold up or lower the bench with minimal effort.

The Pressalit Nursing Bench is perfect for aquatic centres, nursing homes, retirement villages, schools and hospitals.

Features:

- Curved for added comfort and safety
- Manually adjustable back support
- Water collection tray with integral outlet prevents splashing on carer

SPECIFICATIONS

	STANDARD SIZE	JUNIOR SIZE
SKU	3232000	3232001
SWL	200kg	200kg
Overall Length	1910mm	1410mm
Overall Width (from wall)	987mm	987mm
Bed Width	798mm	798mm
Folded Width	290mm	290mm
Height	300-1000mm	300-1000mm
Material	Powder-lacquered stainless steel, aluminium and polyurethane foam	Powder-lacquered stainless steel, aluminium and polyurethane foam
Colour	Grey	Grey
Power Point	10A GPO	10A GPO

Regulation EU 2017/745

