



*Shown in Grey

RELAXER

Dual Motor Lift Chair

Description and Features

- Power lift and recline armchair.
- Chair can be reclined to a trendelenburg position, where the footrest is higher than the head.
- Waterfall style backrest provides superior comfort.
- Padding on each backrest layer can be adjusted.
- Full chaise pad provides additional legrest support and comfort.
- Rise or recline motion stops at any point when button is released.
- Supplied with head cover and magazine pocket.
- Hand held controller is simple to use.
- Dual motor provides independent control of backrest and legrest.
- Available in Fabric or Vinyl and in a variety of colours.
- Made in USA.

Specifications

Product Weight Capacity	170kg
Seat Width	510mm
Seat Depth	535mm
Seat Height	485mm
Back Height	735mm
Arm Height	205mm
Overall Width	930mm
Warranty	12 Months

Fabric Selection

Evergreen Fabric | Dark Brown | Dark Blue/Grey

www.activemobility.com.au
sales@activemobility.com.au
02 9649 2111