

Shoprider 889SL

ASK FOR THE BRAND SHOPRIDER



Available Colours



Active Mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

889SL

www.shoprider.com.au

mobility

The SHOPRIDER™ 889SL is the premium model from the mid size scooters available from the extensive SHOPRIDER range. The 889SL comes with a powerful 4 pole 1.2hp motor, double pillow high back captains seat, adjustable rear suspension and 330mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold down arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

889SL

Length	1300mm
Width	640mm
Height	1200mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.20Hp (4 Pole)
Maximum Speed	10kmph
Battery Capacity (Minimum Amp)	40Amp
Maximum Range #	up to 40km
Wheel Size	330mm
Turning Radius	1440mm
Weight (with battery)	99kg

FEATURES

- Large Captain's Seat
- Double Pillow Comfort Seat
- Adjustable Steering Tiller
- Adjustable Width Arm Rests
- Fold Arm Rests
- Sliding Rotating Seats
- Powerful 4 Pole Motor
- Front Shopping Basket
- Modesty Panel
- Variable Seat height
- Horn
- Battery Level Warning Light
- Front & Rear Bumper Pads
- Front & Rear Suspension
- Front & Rear Lights
- Hazard Lights
- Front Lights
- High / Low Speed switch
- Anti Tips Wheels

Comprehensive Warranty


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



*(Single user)

Unit 8, 110-120 Silverwater Rd
SILVERWATER NSW 2128

Tel: 02 9649 2111

E: Sales@activemobility.com.au

www.activemobility.com.au



Active Mobility

Updated December 2017