

# Shoprider Crossover

ASK FOR THE BRAND SHOPRIDER



**GK10**

Available Colours



Active Mobility

RELIABLE MOBILITY GUARANTEED

# Shoprider

## Crossover

**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au  
*mobility*

The SHOPRIDER™ *Crossover* is a fully portable model that offers exceptional functionality for each user. The Crossover has been designed to offer each user the flexibility of a scooter that can be easily pulled apart for transport or used just as a normal mid size scooter. The Crossover is loaded with features like 265mm wheels, mid back fully rotating captain seat, adjustable steering tiller, multi position lift up armrests & 36amp batteries.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

Length	1190mm
Width	585mm
Height	960mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.2Hp
Maximum Speed	9kmph
Battery Capacity (Minimum Amp)	2 x 36Amp
Maximum Range #	up to 27km
Wheel Size	265mm
Turning Radius	1330mm
Weight (With Battery)	78kg
Weight (Battery Each)	11kg (22kg Total)
Weight (Chair)	16kg
Weight (Rear Transaxel)	17kg
Weight (Main Chassis)	23kg

### FEATURES

Mid Back Captains Seat  
Loop Steering Tiller  
Adjustable Steering Tiller  
Adjustable width arm rests  
Sliding Rotating Seat  
Variable Seat Height  
Battery Power indicator  
Head Lights  
Front & Rear Bump Pads  
Pull Apart Design for Transport  
Shopping Basket  
**Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below


- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

Unit 8, 110-120 Silverwater Rd  
SILVERWATER NSW 2128

Tel: 02 9649 2111

E: Sales@activemobility.com.au

www.activemobility.com.au



Active Mobility

Updated December 2017