

Shoprider

LITTLE BEAUTY

ASK FOR THE BRAND **SHOPRIDER**



Available Colours



Active Mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®
www.shoprider.com.au

mobility

LITTLE BEAUTY

The SHOPRIDER™ **LITTLE BEAUTY** is a versatile portable scooter that offers the performance and features that you expect from a SHOPRIDER™ scooter. This portable scooter is easily disassembled for transport with NO requirement for tools. The **LITTLE BEAUTY** has been designed to offer a number of improved features like LED battery power display, larger leg room, upgradeable battery box for larger batteries and comfortable swivel seat. The reliability of this impressive portable scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

LITTLE BEAUTY

FEATURES

Length	1110mm
Width	540mm
Height	925mm
*Safe Climbing Angle	6 Degrees
*Max Climbing Angle	6 Degrees
Maximum Load Single User*	113Kgs
Charger	2 Amp
Motor Size	0.3Hp
Maximum Speed	6Kph
Battery Capacity (Minimum Amp)	12-21Amp
Maximum Range #	up to 23km
Wheel Size	200mm
Turning Radius	1170mm
Weight (with Battery)	45kg (12 amp Batts)
Weight (Battery Each)	5.6kg (11.2kg total)
Weight (Chair)	8.3kg
Weight (Rear Transaxel)	9.9kg
Weight (Main Chassis)	15.6kg

Mid Back Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Rotating Seat
Folding Arm Rests
Variable Seat height
Front Shopping Basket
Horn
LED Battery Level indicator
Upgradeable Battery Size
Anti Tips Wheels
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.

WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style


Unit 8, 110-120 Silverwater Rd
SILVERWATER NSW 2128

Tel: 02 9649 2111

E: Sales@activemobility.com.au

www.activemobility.com.au

LITTLE BEAUTY



Active Mobility

Updated November 2014