

Wall Mounted Shower Stretcher

The shower stretcher enables people with poor sitting balance to be ergonomically washed under running water. It is attractive, user-friendly and easy-to-maintain with an extensive height range. It is safe for both patients and users. The wall-mounted unit's height is adjustable and has a multifunctional design. By positioning this highly versatile shower stretcher to the correct height, transfers can be easily made from other aids such as beds, stretchers, shower seats and wheelchairs.



Specifications

Product specifications	
Product weight (SWL)	150kg
Weight	85kg
Length	1900mm
Material	Steel, electrolytically galvanised
Warranty	12 months

Features

- 1 folding safety barrier
- Adjustable backrest