

### Soflex Height Adjustable Shower Commode

The Soflex may be utilised for shower and toilet transfer of patients, access to toilet or shower needs, body care, complete safety - all this courtesy of a versatile and comfortable commode. The manual or electric elevating function assists carers providing personal care to maintain a safe posture. The tilt-in-space function provides postural support for clients with limited upper body control.



#### Specifications

Product specifications	
Product weight capacity	135kg
Product length	740mm
Product width	615mm
Product weight	35kg
Seat height	500 – 1505mm
Material	Rislan coated steel - PVC
Colour	White / maroon
Warranty	12 months

