

MiniLift125

SystemRoMedic™

Unique solution

*Small and
convenient*

*Activates
the user*



MiniLift125 is a mobile sit-to-stand lift which, together with a specially adapted rising vest, enables users with reduced function to rise up. The construction of MiniLift125 offers a totally unique and very comfortable sit-to-stand experience. It is also small and convenient and very easy to manoeuvre, which, in combination with the possibility to adjust the settings to the user, makes it very useful in most spaces and situations.



Active Mobility

MiniLift125

MiniLift125 is a mobile sit-to-stand lift which has been developed to, as gently and comfortably as possible, assist the user when rising from a sitting to a standing position. The handles and the side supports offer multiple options for the user to place his/her hands. The stable anti-slip foot plate is very low and the lower leg supports as well as the lifting arm, can be adjusted to fit the user. When MiniLift125 is combined with the appropriate lifting accessories the user gets support under the feet, for the front of the lower legs and behind the back, which provides for a safe and secure sit-to-stand procedure. The construction safely moves the user forwards and upwards in a natural movement pattern and, at the same time, the leg muscles and balance is exercised. MiniLift125 is manufactured in steel for stability and strength and is suitable for users weighing up to 125 kg/275 lbs.



TECHNICAL INFORMATION

Material	Powdercoated steel
User length	Min. 140 cm / 56" , max. 200 cm / 80"
Length of base	900 mm / 36"
Width of base, inner measurement	445 - 745 mm / 17,8" - 29,8"
Width of base, outer measurement	590 - 890 mm / 23,6" - 35,6"
Height	980 mm / 39,2"
Height of footplate	70 mm / 2,8"
Castor size	Ø 80 mm
Lifting speed	3,5 cm/s (1,4 "/s) without load
Emergency lowering	Manual and electrical
Weight	39 kg / 85,8 lbs
Safe working load	125 kg / 275 lbs



Productspecification Art. No
MiniLift125 400641434

ACCESSORIES

MiniLift125 is used together with ThoraxSling, a specially adapted sit-to-stand vest, which is available in many different sizes. ThoraxSling is very easy for the carer to fit and it provides support and security for the user when rising and for standing exercises. It also makes toileting chores easier. MiniLift125 and ThoraxSling provides a very gentle and pleasant sit-to-stand experience. The movement of the lift arm is designed so that the pulling motion ends comfortably in the lower back, which prevents the vest from riding up under the arms.

Disposable material for maximum hygiene and safety

ThoraxSling, with seat support is also available in a strong and durable non-woven material which is disposed of after use, when it gets soiled or when the user does not need it anymore, instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacterial.

ThoraxSling

SystemRoMedic™

ThoraxSling	Art. No	Specification
	45500003	ThoraxSling, XS
	45500004	ThoraxSling, S
	45500006	ThoraxSling, M
	45500007	ThoraxSling, L
	45500008	ThoraxSling, XL
	45500009	ThoraxSling, XXL
	SWL 300 kg / 660 lbs	

ThoraxSling, with seat support

ThoraxSling, with seat support	Art. No	Specification
	45600004	ThoraxSling, with seat support, S
	45600006	ThoraxSling, with seat support, M
	45600007	ThoraxSling, with seat support, L
	45600008	ThoraxSling, with seat support, XL
	SWL 300 kg / 660 lbs	

ThoraxSling, disposable

ThoraxSling, disposable	Art. No	Specification
	45690004	ThoraxSling, with seat support, disposable, S
	45690006	ThoraxSling, with seat support, disposable, M
	45690007	ThoraxSling, with seat support, disposable, L
	45690008	ThoraxSling, with seat support, disposable, XL
	5 pcs	
	SWL 200 kg / 440 lbs	

ThoraxSling, disposable	Art. No	Specification
	70200033	CalfStrap, for MiniLift125 with serial numbers from 125001653

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

Transfer, between two locations. **Positioning**, within the same location.

Support, for mobility. **Lifting**, both manual and mechanical.



Active Mobility

Active Mobility Systems
Unit 8, 110-120 Silverwater Rd
SILVERWATER NSW 2128
Tel.: 02 9649 2111
Fax: 02 9649 8506
E: sales@activemobility.com.au
www.activemobility.com.au