

TA front wheel drive wheelchair with unique maneuverability

TA iQ FWD

Vare number: 97070

Technical data:

Wheelchair Type: CLASS

Max speed forward: 10.00 km / h

Max speed backwards: 5 km / h

Max user weight: 140 Kg

Lifting height: 30 cm

Tilt in space: 45 degrees

Driving distance: 35-40 km

(optimal conditions)

Length without leg supports: 82 cm

Width without seat: 63 cm

Seat Mounting height:

38.5 cm (the top plate)

40 cm (for plate bottom)

41 cm (strap bottom)

Turning radius: 65 cm

Tires (rear wheel): 3.00-8"

Tires (drive wheel): 200-50



- TA iQ FWD is a highly manoeuvrable front wheel drive wheelchair that is easy to manage and get around with both outdoor and indoor.

The small drive wheel is mounted on gas shocks with an adjustable coil spring that provide great soft suspension and a quiet ride over bumps.

- TA iQ FWD comes as other wheelchairs in TA iQ product line with a lift height of as much as 30 cm via a discrete lifting column.

TA iQ FWD can be supplied with a color display control box with straightforward icons that makes it easy to scan the menu, or with a control box without display with symbols for all of the chairs functions.

- TA iQ FWD can be programmed specifically for the user's needs, like all models in the TA iQ product line.

TA iQ FWD can be supplied with Varilite seating system, which includes many different backs and pillow types in various sizes. It Can also be fitted with the most popular seating systems. Seat plate and back support bracket can be adjusted from 38 cm to 55 cm width which allows great flexibility and the possibility of redeploying the chair for other users.