



Active Mobility

Theorem Concepts

LATITUDES

RECLINER AND LIFTCHAIR COLLECTION

5 YEARS LIMITED WARRANTY

www.activemobility.com.au
Sales@activemobility.com.au
02 9649 2111

HOXTON

Single motor lift chair

Fabric Selection

-  Stonewash - Glove
-  Stonewash - Mossy



Specifications

Weight Capacity	113.3 kgs (250 lbs)
Seat Width	49.5 cm (19.5")
Seat Depth	50.8 cm (20")
Chair Width	82.5 cm (32.5")
Chair Depth	99.0 cm (39")
Chair Height	101.5 cm (40")
Operation	Handwand

Description and Features

This single motor lift chair with lay-flat recliner function is the contemporary solution to a functional need. Beautifully upholstered, the Hoxton sets a new standard for the single motor lift chair.

- Single motor lift
- Lay-flat recline
- Handwand



* Shown in Stonewash - Glove

SMYTH

Single motor lift chair

Fabric Selection

-  Stonewash - Dune
-  Stonewash - Amber



Specifications

Weight Capacity	113.3 kgs (250 lbs)
Seat Width	58.4 cm (23")
Seat Depth	53.3 cm (21")
Chair Width	80.0 cm (31.5")
Chair Depth	100.3 cm (39.5")
Chair Height	104.0 cm (41")
Operation	Handwand

Description and Features

Beautifully streamlined, the Smyth has a distinctly Scandinavian feel, while still retaining the smooth single motor lift and lay-flat recline functions. The Smyth is elegantly appointed for unrivalled comfort.

- Single motor lift
- Lay-flat recline
- Handwand





* Shown in Stonewash - Dune

QUIN

Lay-flat recliner with headrest

Fabric Selection

-  Mosaic - Truffle
-  Mosaic - Fudge



Specifications

Weight Capacity	136.0 kgs (300lbs)
Seat Width	58.4 cm (23")
Seat Depth	53.3 cm (21")
Chair Width	80.0 cm (31.5")
Chair Depth	105.5 cm (41.5")
Chair Height	104.0 cm (41")
Operation	Handwand

Description and Features

The Quin is a sleek and contemporary lay-flat recliner with an extended footrest for the ultimate relaxation experience. With the added feature of an independent motorised headrest, the Quin allows you to assume the perfect seating position.

- Two motor lay-flat recliner
- Independent headrest function
- Extended footrest
- Handwand with USB port



* Shown in Mosaic - Truffle

LUDLOW

Dual motor lift chair with headrest and lumbar

Fabric Selection

Stonewash - Dove

Stonewash - Fudge



Specifications

Weight Capacity	158.7 kgs (350 lbs)
Seat Width	53.3 cm (21")
Seat Depth	53.3 cm (21")
Chair Width	83.8 cm (33")
Chair Depth	104.2 cm (41")
Chair Height	105.5 cm (41.5")
Operation	Handwand

Description and Features

The Ludlow dual motor lift chair with headrest and lumbar offers you premium comfort with the unique feature of four individual functions, providing ultimate versatility and superior quality.

- Independent headrest, back, leg and lumbar function
- Extended leg rest
- Handwand with USB port



* Shown in Stonewash - Dove

CROSBY

Three motor lay-flat recliner with headrest and lumbar

Fabric Selection

Stonewash - Gun Metal



Specifications

Weight Capacity	136.0 kgs (300 lbs)
Seat Width	53.3 cm (21")
Seat Depth	53.3 cm (21")
Chair Width	85.0 cm (33.5")
Chair Depth	105.5 cm (41.5")
Chair Height	106.7 cm (42")
Operation	Side Button

Description and Features

The ultimate contemporary lay-flat recliner, the Crosby delivers unrivalled comfort with multiple positional choices thanks to its independent headrest and lumbar functions.

- Three motor lay-flat recliner
- Independent headrest and lumbar functions
- Side buttons with USB port



* Shown in Stonewash - Gun Metal

MERCER

Dual motor lift chair with headrest and lumbar

Fabric Selection

Badlands - Saddle

Badlands - Eclipse



Specifications

Weight Capacity	158.7 kgs (350 lbs)
Seat Width	50.8 cm (20")
Seat Depth	54.6 cm (21.5")
Chair Width	86.4 cm (34")
Chair Depth	104.0 cm (41")
Chair Height	106.7 cm (42")
Operation	Handwand

Description and Features

This dual motor lift chair brings comfort to a new level. The larger than life Mercer is perfect for relaxing and taking time out. With the inclusion of an independent headrest and lumbar, you can adjust the Mercer's settings to form the perfect position.

- Dual motor lift and lay-flat recliner
- Four motors with independent headrest and lumbar
- Extended footrest
- Grained leather look
- Handwand with USB port





* Shown in Badlands - Saddle

SHOREDITCH

Lay-flat recliner with headrest, lumbar and Bluetooth

Fabric Selection

-  Picasso - Walnut
-  Picasso - Midnight



Specifications

Weight Capacity	136.0 kgs (300 lbs)
Seat Width	50.8 cm (20")
Seat Depth	54.6 cm (21.5")
Chair Width	86.4 cm (34")
Chair Depth	104.0 cm (41")
Chair Height	106.7 cm (42")
Operation	Side Button + Bluetooth

Description and Features

The ultimate lay-flat recliner with independent headrest and lumbar, the Shoreditch epitomises luxury. Its rich upholstery and styling gives the Shoreditch the look of a classic recliner with a modern edge, thanks in part to the inclusion of Bluetooth.

- Lay-flat recliner with extended footrest
- Three motors with additional headrest and lumbar
- Side buttons with USB port
- Compatible with iOS devices
-  Bluetooth



* Shown in Picasso - Walnut

STUDIO

Dual motor lift chair with headrest and lumbar

Fabric Selection

-  Santa Barbara Latte
-  Santa Barbara Bark



Specifications

Weight Capacity	158.7 kgs (350 lbs)
Seat Width	57.1 cm (22.5")
Seat Depth	54.6 cm (21.5")
Chair Width	96.5 cm (38")
Chair Depth	106.7 cm (42")
Chair Height	106.7 cm (42")
Operation	Handwand

Description and Features

The Studio dual motor lift chair with headrest and lumbar offers you premium comfort with the unique feature of four individual functions, providing ultimate versatility and superior quality.

- Independent head, back, leg and lumbar functions
- Extended leg rest
- Handwand with USB port



* Shown in Santa Barbara - Latte

ETON

'Box in box' dual motor lift chair with headrest

Fabric Selection

-  Stonewash - Amber
-  Stonewash - Dune



Specifications

Weight Capacity	180 kgs (400 lbs)
Seat Width	53.3 cm (21") 52.0 cm (20.5")
Seat Depth	87.6 cm (34.5")
Chair Width	95.3 cm (37.5")
Chair Depth	106.7 cm (42")
Chair Height	Handwand
Operation	

Description and Features

Aside from its stylish appearance, the Eton also boasts the new integral 'box in box' mechanism. This fully enclosed feature means that nothing can inadvertently obstruct the chair in its fully lifted position. The Eton also includes an independent headrest function for enhanced comfort.

- New 'box in box' technology
- Two motor lift recliner
- Independent headrest function
- Handwand with USB port



* Shown in Stonewash - Amber