

LUDLOW

Dual motor lift chair with headrest and lumbar

Stonewash (Dove)

Stonewash (Fudge)



Shown in Stonewash (Dove)

Specifications

Weight Capacity	158.7 kg (350 lbs.)
Seat Width	53.3 cm (20.9")
Seat Depth	53.3 cm (20.9")
Seat-to-Floor	46.9 cm (18.4")
Chair Width	83.8 cm (32.9")
Chair Depth	104.2 cm (41")
Chair Height	105.5 cm (41.5")
Operation	USB Hand Wand

Description

The Ludlow dual motor lift chair, featuring reclining headrest and lumbar, provides premium comfort with the unique inclusion of four individual functions. With the touch of a button, adjust the chair into your desired position to deliver superior support.

Quick Information

- Dual motor lift chair
- Independent headrest, backrest, legrest and lumbar function
- Extended footrest