

MARCOS

Dual motor lift lay-flat recliner with headrest and lumbar

Mosaic (Fudge)

Mosaic (Charcoal)



Shown in Mosaic (Fudge)

Specifications

Weight Capacity	158.7 kg (350 lbs.)
Seat Width	57.1 cm (22.4")
Seat Depth	55.8 cm (21.9")
Seat-to-Floor	49.5 cm (19.4")
Chair Width	88.9 cm (35")
Chair Depth	105.4 cm (41.4")
Chair Height	107.9 cm (42.4")
Operation	USB Hand Wand

Description

The Marcos dual motor lift chair sets a new standard for relaxation, with the inclusion of an independent headrest, legrest and lumbar function, two discreet cup holders and an easy-to-use USB hand wand control.

Quick Information

- Independent head and extra-wide legrest function lumbar function
- Hidden cup holders
- Hand wand with USB connectivity
- Your choice of two distinctive colours